100*60mm 正反面彩色印刷 层叠折

RIKI SCULPTOR

User Manual



1 Instructions

The RIKI SCULPTOR is a versatile electric fascia device that combines heat therapy, EMS, vibration, and a muscle scraper in one unit. It functions as a muscle relaxer, medical-grade scraper, and fascia massager, providing deep-tissue therapy, physiotherapy support, and facial relaxation.

2 Inside the Box



RIKI SCULPTOR | Charging Cable | Instruction Manual

3 Product Description



4 Safety Information

Important: Improper use of this device may result in injury or damage.

1. Use with care to avoid burns or injuries.
2. Always switch off the device immediately after use and before cleaning to prevent damage.
3. This product is not suitable for individuals with physical, sensory, or cognitive disabilities, or those lacking experience or knowledge (including children), unless supervised by a responsible adult.
4. This device is not water-resistant—do not immerse it in water.
5. The device charges via a 5V USB port, with a full charge taking about 1.5 hours.
6. Keep the device away from high temperatures and humid environments to avoid damage.
7. Avoid dropping or hitting the device, as it may cause operational issues.
8. Store the device out of reach of children.
9. If the device malfunctions, stop using it immediately and contact customer support. If you experience discomfort or pain during use, stop immediately and consult a healthcare professional. This product is not recommended for individuals with skin conditions such as dermatitis, eczema, or swelling.

Heads-up:

Do not use on the following persons:

1. Patients with pacemakers or severe heart rhythm disorders
2. Individuals with cognitive impairments or skin sensory paralysis
3. Infants, pregnant women, and women within 60 days postpartum

Do not use on the following body areas:

1. Areas with metal or plastic implants
2. Skin with open ulcers or wounds
3. The heart, both sides of the carotid artery, the front of the head, and the throat

5 Charging

Connect the device to the charging cable and adapter, then plug it into a household power outlet. While charging, the battery icon will flash. Once fully charged, the icon will stop flashing. The device cannot be powered on while charging.

6 Functional Modes - How to Use

Mode 1:Therapeutic Warmth +Vibration

1. Press and hold the Power Button for two(2) seconds to power on the RIKI SCULPTOR.
2. Single press the Power Button to change the intensity levels from 1~3.
3. Long press the M button to power off/on vibration.

Mode 2: Vibration +EMS

1. Press and hold the Power Button for two(2) seconds to power on the RIKI SCULPTOR .
2. Single press the M button to enter Mode 2.
3. Single press the Power Button to change the intensity intensitylevels from 1~3.
4. Long press the M button to power off/on vibration.

Mode 3: Therapeutic Warmth +Vibration +EMS

1. Press and hold the Power Button for two(2) seconds to power on the RIKI SCULPTOR.
2. Single press the M button twice to enter Mode 2.
3. Single press the Power Button to change the intensity levels from 1~3.
4. Long press the M button to power off/on vibration.

7 Gua Sha Body Techniques

Gua Sha Pushing & Scraping:



Back | Leg | Hand

8 Product Specifications

Product Name	RIKI SCULPTOR
Charging Voltage	5V
Charging Current	1A
Charging Method	USB Charging
Product Weight	296g
Materials	ABS+PC+Zinc Alloy+Silicon
Battery Capacity	1200mAh
Waterproof Rating	IPX0
Product Size	6.42 × 2.22 × 1.18in (163*56.5*30mm)

⚠ WARNING: FOR INDOOR USE ONLY

This device complies with Part 15 of the FCC Rules. Operation is subject to thefollowing two conditions:

(1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE:This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency, energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipmentt does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the used is encouraged to try to correct the interference by one or more of the following measures:

-Reorient or relocate the receiving antanna.
-Increase the separation between the equipment and receiver.
-Connect the equipment into an outlet on a circuit different from that towhich the received in connected.
-Consult the dealer or an experienced radio TV technician for help.

CE FC

RIKI LOVES RIKI™ and RIKI LOVES RIKI™ are trademarks of GLAMCOR GLOBAL LLC.

GLAMCOR GLOBAL LLC
99 Wall Street, 3103
New York, NY 10005

www.rikilovesriki.com
Customer Service: realperson@glamcor.com

©2026 All rights reserved.
Made in China