RIKI BABY FACE LED MASK

User Manual




Table of Contents

What's Inside	1
How Different Light Therapy Works	2
Product Specifications	3
Description for Buttons and Indicators	4
How To Use	5
Recommendation of Use	6
Contradictions and Adverse Reactions	7
Precautions and Travel Tips	8
Intended Use and Instruction for Use	9
Light Sensitivity Test and Preparation	10
Operation Instructions	11
Maintenance	12
Trouble Shooting	13

RIKI BABY FACE LED MASK

What's inside?



Product Description

The RIKI Baby Face is a revolutionary LED light therapy mask designed to give you the radiant, youthful complexion you've always desired. Combining the power of four different light modes (red, red and infrared, blue, and purple), this state-of-the-art skincare device offers an array of benefits to revitalize and rejuvenate your skin.

How Does Red, Blue, Infrared+Red and Purple Light Therapy Work?

Red Light Mode

This mode is perfect for stimulating collagen production, reducing fine lines and wrinkles, and improving skin elasticity. Achieve a firmer, more youthful appearance with consistent use.

Red + Infrared Mode

Combining the power of red and infrared light, this mode is designed to enhance blood circulation, support cell regeneration, and reduce inflammation. Heal and repair your skin from within for a revitalized, glowing complexion.

Blue Light Mode

Say goodbye to acne and breakouts! Blue light therapy targets acne-causing bacteria, reduces inflammation, and balances oil production to promote clear, blemish-free skin.

Purple Light Mode

Harness the benefits of both red and blue light with this mode. Perfect for those looking to combat acne while also addressing signs of aging, purple light therapy offers a comprehensive skincare solution.



RIKI Baby Face Led Mask Specifications

Produce name	
Dimensions	292mm*207mm
Weight	322.5g
LED Bulbs	66 triple-core LEDs
LED Wavelengths	465nm, 630nm, 830nm
Output Power Density	465nm: 15mw/cm^2
	630nm: 15mw/cm^2
	830nm: 15mw/cm^2
Operating Temperature	-10°C to 35°C / -4°F to 95°F
Storing/Transport Temperature	20°C to 60°C / -4°F to 14-°F
Operating Humidity	30% to 80% without condensation
Storage/Transport Humidity	30% to 90% without condensation
Charging Supply and Port	5V-2A, USB-C
Charging Cable	USB-C to USB-A, 1.2 meter / 3.937 feet
Battery	3.7V, 2600mAh Lithium Battery




1 Power/Control Button

2 Battery Indicator

2 Battery indicator:
1. While charging
Battery between 1-20% = Indicator light shows orange.

2. While fully charged:
Indicator light shows solid green.

3. While in normal use:
Battery between 80-100% = Indicator light shows green.
Battery between 1-20% = Indicator light shows orange.



- Make sure your RIKI Baby Face LED Mask is fully charged before use. Takes 4 hours for the battery to be fully charged.
- On clean, dry skin, secure the mask to your face using the Velcro straps, ensuring that it is secure and comfortable.
- Turn on the mask by pressing the controller's power button.
- Select your treatment of either red, blue, infrared+red, or purple light therapy.
- Your 20-minute treatment will begin automatically. Keep the mask on you face for the full 20 minute treatment.
- Turn off the mask by pressing the controller's power button.
- Remove the mask from your face and clean before storing.
- Repeat the treatment three times per week, no more than once per day.

Recommendation of Use

Please read the complete instruction manual before using the RIKI Baby Face LED Mask. It will give you a better understanding of how the product works, and you'll be able to enjoy the best treatment results.

Warning:

1. This device is not for use over the eyes and is not for use on eyelids. Do not look into the light-emitting diodes when using this device.
2. Do not attempt to operate this equipment if there is evidence showing that this product is damaged, or you suspect that this product is damaged. Ask for technical support if you find any damage to the product.
3. Only authorized personnel may service this equipment, contact RIKI LOVES RIKI for servicing or repairing.
4. Do not put body weight on the device. The device cannot withstand a human weight or other heavy external pressure. Do not place the device under clothing, blankets, sheets or towels.
5. Please connect the USB-port to the adaptor before plugging the adapter into to the power supply.
6. Please do not place the RIKI Baby Face LED Mask where it may fall or be pulled into water or other liquid. If the device falls into water or liquid, please unplug immediately.
7. Do not use the RIKI Baby Face LED Mask while you are bathing or taking a shower.
8. Keep the device away from heat and hot surfaces.
9. Do not attempt to disassemble the RIKI Baby Face LED Mask, which may cause damage to the device and void the warranty.
10. Do not let children under 18 years of age use the RIKI Baby Face LED Mask without close supervision.

11. If the device overheats, turn it off immediately.
12. Keep the device away from flammable liquids and aerosols.
13. Use of accessories other than those specified or provided by the manufacturer of this equipment could result in damage or harm to the user.
14. We do not recommend use of the device if you have had a sun burn within two weeks.
15. If you have received any facial fillers or botox within the past month please discuss with your healthcare provider before using the mask.
16. If you have had skin cancer, please discuss with your healthcare provider before using the mask.

Contradictions

(Do not use if any of the following are applicable to you)
1. Do not use the device if you have the Fitzpatrick Skin in Type V and Type VI.
2. Do not use this device if you are photosensitive.
3. Do not use on infants or children.
4. Do not use if you have a medical history of seizures that are triggered by light.
5. Please do not use the RIKI Baby Face LED Mask on the thyroid gland in the neck area.
6. Do not use if you have severe acne.
7. This device should not be used by people who are:
   - Under the age of 18
   - Pregnant or breastfeeding
   - Using photosensitive drugs

Adverse Reactions

If you have any adverse reaction from using this device, stop using it immediately.

Note: If you experience any discomfort, such as skin discoloration, skin damage, or any other concerns develop, stop using the device immediately and consult your physician.

Precautions

1. The intended use of the RIKI Baby Face LED Mask is only as described in the manual.
2. Do not look directly at the light source during the treatment.
3. Do not use the RIKI Baby Face LED Mask if you are sensitive to light.
4. Do not use the RIKI Baby Face LED Mask if you are taking drugs that are known to increase sensitivity to light.
5. Do not use the RIKI Baby Face LED Mask if you are taking steroids or NSAIDs or have received a steroid injection in the last 6 months as you may have a photosensitive reaction.
6. Do not use the RIKI Baby Face LED Mask around other devices that emit electro-magnetic energy in an unshielded manner.
7. Do not scratch or damage the light-emitting surface.
8. Do not use the light of this device to shine on the eyes of others.
9. It will void your warranty if you do not operate the RIKI Baby Face LED Mask according to the manual's instructions.
10. Read all instructions carefully before using this device.
11. Inspect the RIKI Baby Face LED Mask before use in case of the damage.
12. Store the RIKI Baby Face LED Mask according to the suggestions listed in this manual.
13. If you have skin conditions such as Rosacea or Eczema, please discuss with your healthcare provider before using the mask.

Travel Tips

The device's lightweight design makes it a prefect travel companion. There are some travel tips for you:

1. If traveling outside the U.S., do not use the device if the electrical outlet is not within the specifications: AC100-240V, 50Hz/60Hz.
2. Confirm with the airline crew whether the device can be carried and/or used on the airplane.

Intended Use

The RIKI Baby Face LED Mask is intended to emit light in the red region of the spectrum, specifically indicated to treat mild to moderate wrinkles and signs of aging. It also emits light in the blue region of the spectrum, specifically indicated to reduce the bacteria that causes acne and blemishes and improves cell regeneration, skin tone and texture.

Instructions for Use

Check for these situations before use:
Please determine your skin type according to the table below before use, as it will provide more effective instructions for your choice. The device should only be used by Fitzpatrick Skin Types I-IV:

The Fitzpatrick Scale



Type I	Type II	Type III	Type IV	Type V	Type VI
Light, Pale White	White, Fair	Medium, White to Olive	Olive, Moderate Brown	Brown, Dark Brown	Black, Very Dark Brown to Black

Note: The "X" means not suitable. Please don't use the device if you have Fitzpatrick Skin in Type V or Type VI.

Light sensitivity test

Follow the instructions bellow to check if you are light-sensitive:
1. First, turn on the mask (see instructions).
2. Then, apply the mask to the facial area, making sure the light is shining on the skin, and wear for 10 minutes.
3. Wait for 24 hours. If you experience pain or discomfort, or redness that lasts for more than 2 hours, you might be light-sensitive. Please consult a doctor for advice.

First Use and Preparation

1. Take out your RIKI Baby Face LED Mask, Charging Adaptor and Controller from the box.
2. Charge the controller by connecting the DC controller into an electrical outlet (AC100/240V, 50Hz/60Hz).
3. After the controller is fully charged, connect the controller to the mask for the treatment.

Operation Instructions

1. How to prepare for the treatment:
   - Before starting the treatment, thoroughly cleanse and dry your skin, remove any cosmetics that contain a chemical or acidic element that may affect light output, and lead to a decreased or ineffective treatment.
   - Connect the mask and the controller with the USB-cable provided.
2. How to wear the mask properly:
   - Put the mask on your face.
   - Secure the Velcro at the back in whatever position is comfortable, keeping the mask on your face. The Velcro can be secured to whatever fit is comfortable and secure.

3. How to start a treatment:
   - Turn on the mask by pressing and holding the controller's power button to turn the device on. When the mask is turned on with the indicator light on and the beep sound given out, the 20-minute treatment will begin automatically. You can select your treatment of either red, blue, red+blue, or purple light therapy.
4. How long should a treatment last:
   - Once the treatment begins, you should keep the mask on your face for 20 minutes for a full treatment.
5. How to end your treatment:
   - You can turn off the mask in advance by pressing and holding the power button on the Controller.
   - The device will automatically shut down when the 20 minute treatment is over.
6. After a treatment:
   - Remove the mask from your face.
   - Clean the mask (see chapter 6 Maintenance).
   - Store your mask back in the package or in a clean area.

Note:
It is recommended to have a 10-minute treatment, once daily, 3 times a week.
Over-treatment is not advised.
Do not use any makeup or cosmetics that contain chemical or acidic components that may affect light output and lead to a decreased or ineffective treatment.
During the treatment, please use the straps to fix the mask to the face properly.

Maintenance

Cleaning

The mask is reusable and needs to be cleaned after each use.
Please clean the mask after each treatment, using the following method:
1. Use a clean, damp cotton cloth (use water, not other cleaning products) to wipe down the enclosure.
2. Do not oversaturate the cloth; you do not want any liquid to enter the interior of the device.
3. Note: Check the device whether some visible contaminants still exist. If there are visible contaminants on the device, repeat the above steps.
4. Use a clean cotton cloth lightly moistened with alcohol (70%-75%) to clean the enclosure. Like before, be careful to not oversaturate the cloth to ensure that no liquid enters the interior of the device.
5. Visually check the cleanliness of the device and make sure there are no dust, impurities, or any other Pollutant. If necessary, repeat the above step.
6. Drying: please dry the device in the air naturally.
7. Note: Never use abrasive cleaning agents, and never immerse the device in water or other cleaning liquids to avoid damage to the device.

Warning:

a) DO NOT immerse the product in liquid.
b) DO NOT dismantle the unit by yourself.
C) DO NOT drop the device or apply external pressure to the unit since this may damage or weaken the LEDs and other components.
d) If the cord or device is damaged or does not operate normally, contact the distributor. This device may only be serviced, repaired, and disassembled through RIKI LOVES RIKI directly. All units must be returned to the manufacturer for repair, including Warranty repair and Out-Of-Warranty repair, please do not disassemble or repair the returned units by yourself.

Troubleshooting

If the device fails to perform as intended, the following notes will help to identify potential problems with the device and its setup.

Problem	Solution
Device cannot turn on	Ensure that the power supply (controller) of the mask is connected well
	Ensure that the controller of the mask has been fully charged
	Ensure that the controller is unplugged from the power socket
The device suddenly shuts off (out of essential performance)	Ensure you are not using other devices that may interfere with the device
	Ensure that the controller of the mask has been fully charged

Disposal of Electrical and Electronic Equipment Waste (E-waste)

As each city has its regulation of disposal of e-waste, please learn about recycling or disposal programs in your local area through searching "E-waste recycling or disposal" online.
You can also contact your city's department of sanitation (or equivalent office) to learn more about local e-recycling policies.

This RIKI product is proudly designed, tested, distributed, and serviced by GLAMCOR.
• RIKI LOVES RIKI™ and RIKIRIKI® are trademarks of GLAMCOR GLOBAL LLC