RikiMIRACLE

LED ACNE PATCH WITH CHARGING CASE

Instruction Manual

Table of Contents

In the Box & Features	1
How the RIKI Miracle Works	2
Before Using RIKI Miracle	3
Charging	4
Battery Indicators	5
RIKI Miracle Tape	6 -7
Using the RIKI Miracle	8 - 9
Cautions and Contraindications	10 -11
Warnings	12-13
Warranty	13-14
Customer Support	15
RIKI Miracle Technical Specifications	16

1 In the Box



RIKI Miracle Features




2 How the RIKI Miracle Works

The RIKI Miracle combines red and blue light to:

1) Target acne-causing bacteria.
2) Reduce inflammation and redness.
3) Balance sebum production.
4) Enhance skin healing and reduces scarring.
5) Prevent future breakouts.

Benefits: Purple light has advantages of both red and blue skin care light therapy, providing a comprehensive solution for various skin concerns.

Incorporate purple light therapy into your skincare routine for effective, non-invasive acne management and clearer, healthier skin.



3 Before Using RIKI Miracle


Here are some getting started tips.

1) Test the tape on your hand before applying it to your face. A small number of people *may* be sensitive to the tape.
2) Ensure that the treatment area is clean and free from any makeup.
3) Select the correct tape for your treatment.
4) Make sure your RIKI Miracle battery is charged.
5) Do not use the device for more than 40 minutes, as temporary marks may occur.

4 Charging

Insert the supplied charging cable into a standard USB charging port or charging cube. Adapters can be purchased from www.rikilovesriki.com.


Our cable can charge the travel case (A) and the RIKI Miracle directly using the included adapter (B).



5 Battery Indicators

RIKI Miracle will provide three six (6) minute treatments when fully charged.


Do not start the treatment if the ON/OFF button is red as this signifies the battery needs to be charged.



6 RIKI Miracle Tape


Use Tape 1

Smaller blemishes in a concentrated area



Use Tape 2

Larger blemishes in a greater area



7 RIKI Miracle Tape

Why Two Different Shapes of Tape?
We have thoughtfully engineered two variations, each with different cut-out dimensions to accommodate various treatment areas.

This design allows you to position the RIKI Miracle seamlessly over the targeted area while maintaining a buffer, preventing direct contact based on the size of the treatment zone.

Even when addressing open acne, blemishes, or zits, our customizable tape options ensure optimal treatment without further irritation.

Explore the examples to see how these versatile tape options can adapt to your unique needs.

About the tape
Exclusively designed for this product, our adhesive is crafted from skin-friendly materials that prioritize your well-being. It has been rigorously tested on individuals with sensitive skin, ensuring both safety and strength. You can trust in its gentle yet secure hold, making it the perfect choice for users with various skin types.

8 Using the RIKI Miracle

1 Ensure that the treatment area is clean and dry.

2 Choose the tape to accommodate the size of your treatment area.

9 Using the RIKI Miracle

10 Cautions and Contradictions

Do not use RIKI Miracle if you have a history of porphyria, solar urticaria, chronic actinic dermatitis, photosensitivity.

Do not use RIKI Miracle over any suspicious or cancerous lesions.

Do not use RIKI Miracle if you are taking Amiodarone.

Do not use RIKI Miracle if you are undergoing any photosensitizing treatment.

Do not use RIKI Miracle over or near a recent steroid injection.

Do not use RIKI Miracle if you suspect that you have had an adverse/allergic reaction.

11 Cautions and Contradictions

Please read all of the instructions before using the RIKI Miracle. Consult your healthcare provider before using the RIKI Miracle if you have been prescribed any medication for acne.

Pigment darkening caused by visible wavelengths, such as blue light, may be observed in darker skin types and is temporary.

Blue light therapy has not been tested with pregnant women. We recommend that you consult with your GP or midwife before considering the use of RIKI Miracle if you are pregnant.

It is not recommended for children under 13 years old. If you are unsure, please consult a healthcare professional.

12 Warnings

The operator is considered the patient.

Do not place RIKI Miracle over the eyes or where the light will shine into the eyes during treatment.

Do not look directly at the LED lights when it is switched on.

Always keep the RIKI Miracle away from water and moisture.

Never insert objects into the openings of the RIKI Miracle.

Never use the RIKI Miracle if it appears damaged or if any of the power cables appear damaged.

Never disassemble the RIKI Miracle.

Do not allow any fluid to enter the RIKI Miracle at any time. Keep the device out of reach of children. The device contains a cable that may pose a risk of strangulation.
The RIKI Miracle Strips should always be carefully removed after use to help avoid damaging the skin.

13 Warnings

Cleaning
The cleaning of the RIKI Miracle device should be performed by the operator in a home healthcare environment. Waste Electrical and Electronic Equipment (WEEE): Do not dispose of the device with normal domestic waste at the end of its life. Please take it to an official recycling point for proper disposal. To recycle this device, you can return it to RIKI LOVES RIKI.

Note: The RIKI Miracle is not waterproof.
If the RIKI Miracle device requires cleaning, use a soft dry cloth or a cleansing wipe. Do not allow any fluid, including cleaning detergent, to enter any part of the RIKI Miracle device. If this happens, do not use the device.

For hygiene reasons, do not share your RIKI Miracle device with anyone else. The cleaning of the RIKI Miracle device should be performed by the operator in a home healthcare environment.

14 Warranty

RIKI Miracle treatment for acne is warranted to be free from defects in material and workmanship for a period of 12 months from its original purchase date under normal use. This warranty applies only to the original retail purchaser with valid proof of purchase and when purchased from an authorized retailer or reseller.

Please note that this warranty is valid only when the RIKI Miracle device is used in conjunction with authorized RIKI Miracle strips.

If the product becomes defective within the warranty period, please contact the retailer from whom you purchased the device.

This warranty does not cover products damaged by:
- Accidents, misuse, abuse or alteration.
- Use with unauthorized accessories or use other than as instructed.
- Connecting to incorrect current and voltage.

15 Customer Support

We understand the importance of your Riki Loves Riki product, and we're here to ensure you have the best possible experience with it.

Should you encounter any issue, we're just an email away. Reach out to us at help@rikilovesriki.com and our dedicated customer solutions team will swiftly assist you with your warranty claim.

For more comprehensive information and additional support, we welcome you to visit our Customer Solutions page at:

http://www.rikilovesriki.com/pages/customer-solutions.

Our goal is to promptly resolve your concerns and ensure your complete satisfaction.

Remember, at RIKI LOVES RIKI, we are committed to delivering quality products and exceptional service every time!

16 Product Specifications

RIKI Miracle Specifications

Product Name	RIKI Miracle
Dimensions	65mm*36mm*26mm
Weight	322.5g
LED Bulbs	6 double-core LEDs
LED Wavelengths	415nm, 630nm
Output Power Density	215nm: 20mw/cm² 630nm: 20mw/cm2
Operating/Storing Temperature	0°C to 35°C /32°F to 104°F
Operating Humidity	30% to 80% without condensation
Storage/Transport Humidity	30% to 90% without condensation
Charging Supply and Port	5V-2A, USB-C
Charging Cable	USB-C to USB-A
Charging Cable Length	1 meter / 3.937 feet

RIKI LOVES RIKI™ and RIKI❤RIKI™ are trademarks of GLAMCOR GLOBAL LLC

This RIKI product is proudly designed, tested, distributed, and serviced by GLAMCOR

US and other foreign patents pending.

©2023 All rights reserved.
Made in China